



PROFESSIONAL SYSTEM SEGREGATION AND STORAGE DRUGS



✉ PL 62-070 Dąbrówka k/Poznań, ul. Widok 11

☎ (61) 816 13 11 📠 (61) 415 16 75

🌐 biuro@techniapteka.pl

💻 www.techniapteka.pl

GRUPAMM

Thanks to deep stacking drawer systems STEEL BOX, the time necessary to localize a given medicine is shortened to a minimum, which allows for quicker service and devoting more time to a patient. Systems of organization and storage of medicinal drugs STEEL BOX allow for maximal decreasing of storage areas. Casing well-fitting the existing furniture in the pharmacy or client's expectations. Quality of workmanship confirmed by a 5-years warranty.

All metal elements were covered with the best quality powder paints. Antistatic surfaces prevent dust gathering.

Stable construction and big storage capacity. Light and exceptionally rigid sides with profiled notches for compartments.

Occupational safety and hygiene. Metal grid at the bottom of drawers prevents dust from gathering. Automatically closing drawers minimize the risk of suffering an injury.

Quicker access to articles, thanks to better designed storage area.



Good organization. Easy identification of products, thanks to unique compartments and marking system, which minimizes the risk of making a mistake while dispensing medicines.

Silent work of driving mechanisms.


More room. In view of storage capacity two STEEL BOX columns correspond to six meters of traditional stack columns.

DIMENSIONS OF STEEL BOX SYSTEM [mm]

	A	B	C
single	1150	121	2×172
double	1150	245	2×172



COMBINATION OF DRAWERS USED THE MOST OFTEN



m	28	26	24	16
m ²	4,82	4,47	4,13	2,75
kg	152,4	145,8	139,2	112,8

